

DiluFlow®

Gravimetric dilutors for solid samples



The best solution for your dilutions!

ambifood®




Our quality for your lab

- Full range of products for microbiological analyses
- Leadership for innovative products
- Worldwide distribution network
- Designed and Made in France




BagSystem®

for the best sample preparation



Why use a gravimetric dilutor?

DiluFlow® gravimetric dilutors enable the automatic dilution of the sample with perfect reliability. See the following advantages:



Time-consuming ■

Inaccurate: ■
Direct impact
on final result

Not reproducible ■

Not secure: ■
The bag is not maintained

Not automatically traceable ■

■ **Quick:** 8 seconds* for 225 mL

■ **Accurate:** < 1 % of sample weighing
As little as < 1 % of dispensing
(in standard or accurate mode)

■ **Reproducible**

■ **Secured:** The bag is maintained
with BagOpen® bag opener

■ **Traceable:** with LIMS, printer or Excel™

* Dilution time for a 25 g sample diluted at 1/10 on a DiluFlow® Elite with booster kit in fast mode

Easy-to-use efficient & connected

3 models for optimal dilutions



	DiluFlow®	DiluFlow® Pro	DiluFlow® Elite
Ultra low profile	✓	✓	✓
DripTray	✓	✓	✓
Robotic arm	-	✓	✓
Jumbo kit option	-	✓	✓ ^(*)
Programmable	-	-	✓
Full traceability	-	-	✓
Up to 6 pumps	-	-	✓
Multi-Dispensing	-	-	✓

Superior quality

Watson Marlow™ pumps



Traceability

Printer, USB stick⁽¹⁾, Excel™, Open Office™, Libre Office Calc™, CSV, LIMS⁽¹⁾



Up to 6 diluents

With 4 extra connectable pumps⁽¹⁾

IN COMPLIANCE WITH

ISO
7218

ISO
6887-1

FDA BAM
Bacteriological Analytical
Manual

304 L stainless steel casing

For easy cleaning and long-lasting

Robotic arm⁽²⁾

Easy sample dropping
Avoids contamination
Saves time

Demo video
Scan me !



GeckoGrip

No more adhesive gums!
Ultra adhesive, hygienic and long lasting

PATENTED

Ultra low profile

Perfect for a laminar flow
Ergonomic: no more
musculoskeletal disorders

Kit booster

2 times faster

Only 30 cm!

Removable DripTray


Easy to keep your DiluFlow® clean
PATENTED

High accuracy weighing scale

Weighing accuracy up to 99 %

(1) Only on DiluFlow® Elite (2) Only on DiluFlow® Pro and DiluFlow® Elite (3) Only on DiluFlow® Elite 5 kg

Unique features




World's fastest dilutor:

With a booster kit, the Dilu**Flow**® are the fastest dilutors on the market. Only 8 s⁽³⁾ to dispense 225 mL!

Robotic arm⁽²⁾ & nozzle protection plate

The robotic arm enables you to drop the sample safely and easily without contamination risks. Tested more than 1.3 million times, its reliability is guaranteed.

The protective plate protects the nozzle from any contact with the sample.



Ultra low profile

The Dilu**Flow**® have the lowest bag opening on the market, less than 30 cm!

This allows a better comfort and ergonomics, especially when you use the DiluFlow in a biosafety cabinet.

GeckoGrip gripping system

Inspired by a gecko's paw, this new GeckoGrip material will hold the bag without pressing hard. GeckoGrip sticking properties apply in dry and wet surfaces.

GeckoGrip long lasting material can be washed with soap and alcohol.

PATENTED





LightCode

An intuitive indicator of the status of the DiluFlow®. The degree of accuracy of the finished dilution is displayed: green > 98 %, orange > 95 % or red < 95 %. It will also turn blue if the scale is not stable, thus telling you to wait for stabilization.



Green:
Accuracy > 98 %



Orange:
Accuracy > 95 %



Red:
Accuracy < 95 %



Blue:
Scale stabilizing

DripTray

In case culture media is spilt on the unit, the new patented DripTray can be removed to be rinsed, avoiding a tedious cleaning of the unit.

PATENTED



Multi-Dispensing⁽¹⁾


This unique patented mode allows you to serially dispense in small containers such as 9 mL tubes, without having to weigh each tube, by doing just one calibration on the scale.

PATENTED

Magnetic BagOpen®

The 2 points magnetic attachment of the **BagOpen®** guarantee a stable platform for better measuring accuracy, and provide a natural protection of the weight sensor when removing the bag.


PATENTED SHAPE



(1) Only on DiluFlow® Elite (2) Only on DiluFlow® Pro and DiluFlow® Elite (3) Dilution time for a 25 g sample diluted at 1/10 on a DiluFlow® Elite with booster kit

Multiple ways to use a DiluFlow®

☞ The lab operators like the low profile of the DiluFlow® platform ☞



New standards make it necessary to work under a laminar flow for the dilution work, and this can create operational problems.

Our lab operators like the ultra-low profile of the DiluFlow® platform because they do not need to bend themselves anymore to insert the sample in the bag.

My DiluFlow® can handle up to 5 000 g load if you add the Jumbo kit option, to dilute in large bags. It is useful because you do not need to buy an extra machine for large samples.

■ Brigitte, lab technician

☞ The DiluFlow® Elite, in the 1 kg version, has an accuracy of 0.01 g ☞



In the cosmetics industry, the raw materials cost a fortune, sometimes more than 10 000 euros per kilo ! That is why we only work with small amounts, usually 1 g.

The DiluFlow® Elite, in 1 kg version, has an accuracy of 0.01 g, which enables us to get $\pm 2\%$ weighing accuracy. Moreover we use the tongs platform option which are efficient pliers to hold our tubes, to directly place the tube under the liquid flow every time.

■ Dylan, lab manager

☞ Accurate, user-friendly and reliable... ☞



We analyze between 10 and 20 samples per day in our lab.


We looked at all dilutors on the market and most of them looked far too complex to use or too sophisticated.

We are now using a DiluFlow® which matches our requirements completely. The keyboard is really easy-to-use; only 2 keys to press for a dilution.

■ Alain, microbiological coordinator



☞ 3-in-1 equipment ☞




What I really like about my DiluFlow® Elite is the Multi-Dispensing function. I use it to fill 9 mL dilution tubes. Before that, I was using a peristaltic pump, but each time, I had to calibrate the pump with a weighing platform to get $\pm 2\%$ accuracy.










With the DiluFlow® Elite, you just put the first tube on the platform and the calibration is automatic. You can then directly distribute the doses in the tubes with the HandyGun, without having to weigh each tube. It is as if I had 3-in-1 equipment: a peristaltic pump, a scale and a dilutor!



■ Dominique, quality control engineer

Traceability



	DiluFlow® DiluFlow® Pro	DiluFlow® Elite
Available ports	   Input for pedal RS232 for printer USB for PC connection	    Ethernet Jack USB B USB A
Data export to a printer	 Sample number: 1227E-178-PA Operator: Jean Durant 02/12/2014 10:06:37 Pump: Right Theoretical: 290 g Actual weight: 290,61 Precision: 99,8 % * Printed via the supervision software	 Numéro échantillon : 847610931 24/04/15 15:16:05 Dernier étalonnage : 24/04/15 Opérateur: user 0 Appareil: 200453 Programme: 25g 1/10 dilution Lot: kf/1684 Taux de dilution : 1/10 Masse de l'échantillon : 25.32 g Masse réelle : 255.38 g Précision : 99.7 % * Printed directly from the DiluFlow® Elite
Import and export of computer & network data	Traceability : <ul style="list-style-type: none"> • USB connection to computer and data export via the DiluFlow® supervision software • Excel™ file export from the supervision software 	Full traceability : <ul style="list-style-type: none"> • Data export to Excel™ • Data export to USB key (CSV) • Bi-directional connectivity to the LIMS by Ethernet

How to order



DiluFlow®

Single pump (3 kg)
Double pump (3 kg)

Ref 501 103
Ref 501 203



DiluFlow® Pro

Single pump (5 kg)
Double pump (5 kg)

Ref 502 105
Ref 502 205



DiluFlow® Elite

Single pump (1 kg)
Double pump (1 kg)
Single pump (5 kg)
Double pump (5 kg)

Ref 503 101
Ref 503 201
Ref 503 105
Ref 503 205

DiluFlow® is delivered with 1 dispensing arm, 1 DripTray 400, 1 or 2 pumps, 1 dispensing assembly (GL45 - Ø 6,4 mm) per pump, USB cable, monitoring software, power cord, 1 BagOpen® 400 for dilutors, user's manual, quick user guide, 10 GeckoGrip, free sample bags.
DiluFlow® Pro is delivered with 1 or 2 pumps, 1 DripTray 400, 1 protection plate for tubing nozzle, 1 dispensing assembly (GL45 - Ø 6,4 mm) per pump, USB cable, monitoring software, power cord, 1 BagOpen® 400 for dilutors, user's manual, quick user guide, 10 GeckoGrip, free sample bags.
DiluFlow® Elite is delivered with 1 or 2 pumps, 1 DripTray 400, 1 protection plate for tubing nozzle, 1 dispensing assembly (GL45 - Ø 6,4 mm (Elite 1 kg: Ø 3,2 mm)) per pump, USB stick, power cord, 1 BagOpen® 400 for dilutors, user's manual, quick user guide, 10 GeckoGrip, free sample bags.

Accessories



Printer

Prints all data
Paper 521 110
Labels 521 012
Film 410 130
Ref 410 150



Foot pedal

For hands-free dispensing. Connection with a 3.5 mm Jack.
Ref 507 008



Bar-code reader (1)

For total traceability of the samples.
Ref 522 000



Handy gun (1)

For high precision work and comfortable tube handling.
Ref 561 001



Booster kit

Accelerate the dispensing speed up to 225 mL/8 s.
Ref 511 000



Pump module (1)

External module with single pump + silicone tubings for pump head (inner Ø 6.4 mm).
Ref 505 007



Jumbo kit (2)

BagOpen® 3500
+ 3500 mL DripTray
+ anti-splash nozzle
Ref 505 006



5L bottle (GL45)

5 L autoclavable bottle
Ref 513 001



Independent stand for broth bags

To easily store diluent bags.
Ref 513 015



Integrable stand for broth bags

Integrated to the unit, to hold one broth bag
Ref 505 010




Trocars

316 L stainless steel trocars for broth bags.
Ref 513 020




6 nozzles holder (1)

+ Protection nozzles
Can hold up to 6 nozzles
Ref 505 008




Adjustable distribution system (3)

Ideal to fill various containers
Ref 505 004




Adjustable distribution system (4)

Ideal to fill various containers
Ref 505 005



Weighing platform

To weigh and support a container
Ref 505 002



Tongs platform

To hold a tube (Ø 10 to 60 mm)
Ref 505 003



GeckoGrip

Set of 10 adhesives for BagOpen®
Ref 505 009



PTFE filters 0.2 µm (Ø 25 mm)

Set of 5 autoclavable filters for GL45 bottles
Ref 513 023



Nozzle

Set of 5 stainless steel nozzles for tubings (inner Ø 4.8 to 8 mm)
Ref 513 024



Nozzle

Set of 5 stainless steel nozzles for tubings (inner Ø 1.6 to 3.2 mm)
Ref 513 040



Tube weight (Ø 6 mm)


Set of 5 tube weights
Ref 513 039



Certified standard weight (500 g)


500 g COFRAC calibrated weight, type F1
Ref 513 005

Tubings



Distribution set for single pump **Sales unit**

Ref 513 045	Dispensing assembly inner Ø 3.2 mm	1
Ref 513 043	Dispensing assembly inner Ø 6.4 mm	1




Distribution set for booster kit **Sales unit**

Ref 513 008	Booster kit dispensing assembly inner Ø 3.2 mm	1
Ref 513 041	Booster kit dispensing assembly inner Ø 6.4 mm	1




Silicone tubings for single pump **Sales unit**

Ref 513 047	Silicone tubings inner Ø 3.2 mm	5
Ref 513 044	Silicone tubings inner Ø 6.4 mm	5



Silicone tubings for pump head **Sales unit**

Ref 561 302	Silicone tubings for pump inner Ø 3.2 mm	6
Ref 561 304	Silicone tubings for pump inner Ø 6.4 mm	6



Silicone tubings for booster kit **Sales unit**

Ref 513 009	Silicone tubings for booster kit inner Ø 3.2 mm	5
Ref 513 042	Silicone tubings for booster kit inner Ø 6.4 mm	5

Complementary products



BagMixer®
Lab blenders

Models (mL):
100/400/3500



BagFilter®
Sterile bag with lateral filter

Models (mL):
400/3500



BagPage®
Sterile bag with full-page filter

Models (mL):
100/400/3500



BagLight®
Sterile plain bag

Models (mL):
100/400/3500



BagOpen®
Bag opener/holder

Models (mL):
80/100 - 400 - 3500



BagRack®
Storage rack

Models (mL):
80/100 - 400
3500 (2 or 4 compartments)



BagClip®
Closing clips

Models (mL):
80/100 - 400 - 3500



BagPipet® & BagTips®
Fixed volume pipet & sterile tips

BagTips® Models:
Regular (19 cm)
Jumbo (24 cm)



BagSeal®
Thermal sealing unit

For bags
up to 400 mL



PetriPile®
Petri dish storage rack

Models (ø):
55 - 65 - 90 - 150 mm

Technical specifications

	DiluFlow® 3 kg		DiluFlow® Pro 5 kg		DiluFlow® Elite 5 kg		DiluFlow® Elite 1 kg	
	Single pump	Double pump	Single pump	Double pump	Single pump	Double pump	Single pump	Double pump
Reference	501 103	501 203	502 105	502 205	503 105	503 205	503 101	503 201
BagOpen® height < 30 cm								
Patented removable DripTray								
BagOpen® with GeckoGrip								
Magnetic stand for BagOpen® (patented shape)								
LightCode: color status LED								
Compatible with Booster kit								
304 L stainless steel body								
Stainless casting electromechanical core								
Robotic arm	-							
Compatible with Jumbo kit	-							-
Protection for tubing nozzle	-							
Up to 6 pumps	-			-				
Patented Multi-Dispensing mode	-			-				
32 editable programs	-			-				
Bidirectional connectivity	-			-				
Dilution 1 g at 1/10	-			-		-		
Resolution up to 200 g					0.01 g			
Weighing range	3 000 g				5 000 g		1 000 g	
Maximum weight without support	4 500 g				7 000 g		1 500 g	
Weighing accuracy	0 to 100 g: ± 0,05 g 100 to 3 000 g: ± 0,1%				0 to 100 g: ± 0,05 g 100 to 5 000 g: ± 0,1%		0 to 20 g: ± 0,01 g 20 to 1 000 g: ± 0,1%	
Dilution factor	1/2 to 1/99				2 to 1000, editable fractions			
Dilution accuracy	> 98 %				> 99 %*			
Dilution time for a 25 g sample diluted at 1/10	< 9 s**				< 8 s**			
Dilution time for a 375 g sample diluted at 1/10	-		< 80 s**		< 75 s**		-	
Distribution time for 90 mL	< 7 s**				< 6 s**			
Distribution time for 225 mL	< 10 s**				< 9 s**			
Max number of pumps	2 integrated				2 integrated + 4 external			
Dilution/distribution mode	Fast - Standard - Accurate							
Traceability	Printer / Monitoring software (data export to Excel™, Open Office™, Libre Office Calc™, CSV, printer)				Printer / Export to Excel™ / CSV file export on USB stick / LIMS			
Exported data	Sample weight / Total weight / Dilution factor / Dilution accuracy / Right or left pump + with monitoring software: Sample number / Name of operator / Date / Time				Sample weight / Total weight / Dilution factor / Dilution accuracy / Diluent / Sample number / Diluent batch number / Date / Time / Last calibration date / Name of operator / Name of program			
Connection	RS232 output for printer / USB output for computer / Jack output for foot pedal				2 USB outputs (type A) to connect printer, keyboard, bar-code reader, USB stick / 1 USB port (type B) to connect a computer / Jack output to connect a food pedal / Handy gun / Ethernet connection for LIMS			
Language	English				English, French			
Weight	9,6 kg	12 kg	11,8 kg	14,2 kg	12,6 kg	14,8 kg	12,6 kg	14,8 kg
Dimensions (w x d x h)	31 x 44 x 38 cm	37 x 44 x 38 cm	31 x 42 x 38 cm	37 x 48 x 38 cm	31 x 42 x 38 cm	37 x 42 x 38 cm	31 x 42 x 38 cm	37 x 42 x 38 cm
Patents	GeckoGrip system, Removable DripTray, Multi-Dispensing, Magnetic BagOpen® (patented shape)							
In compliance with	ISO 7218 - ISO 6887-1 - FDA BAM (Bacteriological Analytical manual)							
Warranty	1 year							

* with booster kit in standard or accurate mode **with booster kit in fast mode

Certified production:



ambifood®

Rua Dominguez Alvarez, 44, 4.16 . 4150-801 Porto . Portugal

+ 351 229 962 069 . geral@ambifood.com . www.ambifood.com