




SEROCONVERSION OF DAY OLD COMMERCIAL BROILERS VACCINATED WITH VECTORMUNE-ND*



The **ID Screen® Newcastle Disease Indirect** ELISA is a rapid and reliable alternative to HI for use with recombinant vector vaccines.

The test offers:

- Excellent agreement with HI
- Improved sensitivity
- **Earlier detection** of vaccine antibodies

The IDvet ELISA gave similar results to HI, with a slight increase in antibody levels between D21 and D28, and a sharp increase between D28 and D42. At D35, the IDvet ELISA and Kit A showed 100% and 40% positivity, respectively.

*Study presented by the Scientific Direction of Ceva Animal Health at the 2015 WVPA meeting in Cape Town, South Africa