



PORCINE EPIDEMIC DIARRHEA VIRUS

MONITOR VACCINATION, IDENTIFY AND CONTROL NEW OUTBREAKS

The ID Screen® PEDV Indirect ELISA is a rapid and easy-to-use test for detection of IgG antibodies directed against the PEDV in serum and plasma. Use this kit to:

- Identify and control outbreaks
- Determine antibody prevalence following natural exposure or feedback / vaccination
- Demonstrate the absence of exposure at herd level

- + High performance
- + Easy-to-use: no sample pre-dilution, with colored and ready-to-use dilution buffers
- + Rapid: results in only 90 minutes
- + Data analysis with the IDSoft™ program



Specifications

METHOD	Indirect ELISA
SAMPLE TYPES	Serum and plasma
SPECIES	Swine
PRODUCT CODE	PEDVS-5P
FORMAT	5 plates, 480 tests