



EQUINE INFECTIOUS ANEMIA

AGAR GEL IMMUNODIFFUSION (AGID) TEST FOR THE DETECTION OF EIA ANTI-P26 ANTIBODIES



- Clear, bright lines for easy reading
- OIE-prescribed protocol
- May be purchased with or without agar gel (the IDvet agar gel allows for clearer precipitation lines)
- 24-month shelf-life
- Generates a clean inflexion with OIE reference serum



WITH YOU AT EVERY STEP

ID.vet

Specifications



METHOD	AGID (Coggins-test)
SPECIES	Equine
SPECIMENS	Serum and plasma
ANTIGEN	EIAV p26

Ordering information

PRODUCT CODE	EIA-AGID	EIA-AGID-NOGEL
KIT FORMAT	Includes gel	Gel not included
REACTIONS	200	200

Associated reference serum

PRODUCT CODE	MRI-EIA-AGID
FORMAT	0.5 ml vial (freeze-dried)
DESCRIPTION	Freeze-dried EIA positive horse serum to be used as internal reference material for quality control