

# Anabac<sup>®</sup>

Autoclave deodorant



ambif<sup>∞</sup>d<sup>®</sup>

# Give some air to your lab!

Laboratories are very concerned about hygiene, cleanliness and comfort of work. Try **Anabac**<sup>®</sup>, autoclave deodorant, for a perfumed touch in your lab. **Anabac**<sup>®</sup> autoclave deodorant bring the discrete final touch, the refinement of cleanliness.

## Classic

---



Well perfumed  
All applications  
Ref 320 100

- Mint and eucalyptus fragrance



## Poma

---



Sweet perfume  
Food industry  
Ref 320 200

- Apple fragrance



## Citrus

---



Sweet perfume  
Food industry  
Ref 320 300

- Lemon fragrance



## Floral

---



Well perfumed  
All applications  
Ref 320 400

- Flower fragrance



## Peach

---




Well perfumed  
All applications  
Ref 320 500

- Peach fragrance



1.



*During sample sterilization, bad smell can emerge from the autoclave.*

## Advantages

- Work comfort
- Compatible with any autoclave
- No effect on the autoclaving
- Pleasant perfume
- 100% biodegradable


## Applications

- For any kind of autoclave (max: 200 L), recommended for food industry.

## Specifications


- Pot of 100 capsules
- Pot weight: 160 g
- 5 fragrances:
  - Classic (ref: 320 100)
  - Poma (ref: 320 200)
  - Citrus (ref: 320 300)
  - Floral (ref: 320 400)
  - Peach (ref: 320 500)

2.



*Place one capsule (or 2 if the odor is strong) in the water at the bottom of the autoclave.*

3.



*After a sterilization cycle, the capsule melts, allowing the perfume to be released.*



## Discover also our disinfection range

### Applications

- Laboratory: lab equipment, glassware, benchtops, small instruments...
- Food production sites



### Clinet® 101

*Disinfecting wipes*

Ideal for laboratories, hospitals, samplers, home use.

Ref 310 101



### Clinet® 202

*Disinfecting wipes in individual pack*

Ideal for laboratories, hospitals, samplers, transports, hand hygiene.

Ref 310 202



### Clinet® 303

*Bactericidal and fungicidal spray ready to use*

Ideal for large surfaces, hard to reach places, angle of walls, rough edges.

Ref 310 303



### Clinet® 404 SaniFood®

*Disinfecting wipes for food environment*

For food professionals: disinfection without rinsing of surfaces and materials in contact with food.

Ref 310 404



### Clinet® Anti-algae

*For laboratory water baths.*

1 liter bottle with incorporated doser, Clinet® Anti-algae is used diluted at 0.1% to inhibit the growth of algae in water baths.

Ref 310 400

**ambifood®**

Rua Dominguez Alvarez, 44, 4.16  
4150-801 Porto Portugal

TEL: +351 229 962 069 FAX: +351 229 940 082

geral@ambifood.com

[www.ambifood.com](http://www.ambifood.com)