



CONTROL AVIAN LEUKOSIS IN YOUR FLOCK

The ID Screen® ALV Antigen Capture is a sandwich ELISA for the detection of p27 antigen of the Avian Leukosis Virus in chicken light albumen or cloacal swabs.

ID SCREEN® ALV ANTIGEN CAPTURE

- ➕ Offers the highest analytical sensitivity
- ➕ Detects all Avian Leukosis subgroups, including the new K subgroup circulating in China.
- ➕ High reproducibility, repeatability, stability and robustness
- ➕ Data analysis possible with the IDSoft™ program



Specifications

METHOD	Antigen Capture ELISA based on anti-P27 MAbs	
SAMPLE TYPES	Light albumen or cloacal swabs	
SPECIES	Chickens	
PRODUCT CODES	ALVAG-5P	ALVAG-10P
FORMATS	5 plates, 480 tests	10 plates, 960 tests
PLATE FORMATS	1 x 96-well whole plates	1 x 96-well whole plates